

Patient Electronic Access: View, Download or Transmit

MEASURE DESCRIPTION: At least one patient seen by the MIPS eligible clinician during the performance period (or patient-authorized representative) views, downloads or transmits their health information to a third party during the performance period.

DENOMINATOR: Number of unique patients seen by the MIPS eligible clinician during the performance period.


NUMERATOR: The number of unique patients (or their authorized representatives) in the denominator who have viewed online, downloaded, or transmitted to a third party the patient's health information during the performance period.

Threshold: At least one patient


How to perform View, Download or Transmit using iSmart EHR.

View, Download or Transmit action should be done by patient/patient representative who has access to iSmart Patient Portal. To get information how to provide access to the Portal, please refer to Provide Patient Access measure.

1. A patient is given his/her personal credentials for iSmart Patient Portal. The patient follows the link <https://phr.wchsbc.com/> and enters to the Portal using login and password.



Sign in

 Patient Panel
Smart
EHR

WARNING! UNAUTHORIZED ACCESS PROHIBITED. You must have authorized permission to access this area.

Login or Email:

Password:

Login

[Forgot Password](#)

2. The patient should switch to 'Patient Summaries' page and hit 'New' button.

APPOINTMENTS AMENDMENTS MESSAGES PRACTICE INFORMATION **PATIENT SUMMARIES** ACTIVITY LOG NETWORKING

Gender: F
Date of Birth: 4/17/1951
Race:

Ethnicity:
Preferred Language:
Smoking Status: Ex-moderate cigarette smoker (10-19/day) (finding)

Pharmacy Care
Social Security

Patient

Generation Date	Date From	Date To	Note
-----------------	-----------	---------	------

New

3. Choosing the date for the Clinical Summary and clicking '**Generate**' button, the patient will see newly generated summary.

Current Date: 11/1/2017 8:29 PM

You may filter data by Date range or specific date:

Single date Date range

Select date

Note:

Cancel **Generate**

4. Patient needs to click **'view'** icon to see the generated summary, **'download'** icon to download the summary. In order to transmit summary patient has to click **'transmit'** icon, enter email address to send the summary to and click **'Transmit'** button.

